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<Flex-Beads Preparation>

1. Flex-beads sonicate for 30", vortex, then all transfer to 1.5ml eppendorf.
 2. 10,000rpm, 4min.
 3. Remove the supernatant & add 200ul 1.5XTMAC buffer, then vortex.
 4. Repeat step2), 3) twice (totally wash three times)
(Each residue the supernatant 25ul).
 5. Add 1.5XTMAC buffer 75ul(total 100ul).
- >Final Flex-beads concentration: **2.5*10⁶/ml**, storage at 4⁰C.

<PCR-reaction >

1. After virus RNA extraction, RNA reverse transcribed to cDNA by RT reaction.
2. Then add reagents & cDNA/Plasmid by the table:
- 3 After nested PCR , check the product by gel electrophoresis.

<Hybridization>

1. Add PCR product 10ul, Flex-beads 2ul, and 1.5X TMAC buffer 18ul(each reaction), mix thoroughly.
2. By PCR machine: 96°C 5'(denature)
→37°C 30'(hybridization)
→37°C stand
3. 2,000rpm, 3min.
4. Remove the supernatant & add 75ul 1XTMAC buffer, mix.
5. 2,000rpm, 3min, remove the supernatant.
6. Add 75ul 1XTMAC with 5ug/ml streptavidin-R-phycoerythrin by the table:
7. React at 37°C for 15 min.
8. Mix, then analysis by Luminex 100IS machine.

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1. Whitehead GS, Walker JK, Berman KG, Foster WM, Schwartz DA, Allergen-induced airway disease is mouse strain dependent. *Am J Physiol Lung Cell Mol Physiol*, **285** (2003): L32-L4.
2. Abreu I, Laroche P, Bastos A, Issert V, Cruz M, Nero P, Fonseca JE, Branco J, Machado Caetano JA, Multiplexed Immunoassay for Detection of Rheumatoid Factors by FIDISTM Technology. *Ann N Y Acad Sci*, **1050** (2005):357-63.
3. Gires O, Munz M, Schaffrik M, Kieu C, Rauch J, Ahlemann M, Eberle D, Mack B, Wollenberg B, Lang S, Hofmann T, Hammerschmidt W, Zeidler R, Profile identification of disease-associated humoral antigens using AMIDA, a novel proteomics-based technology. *CELLULAR AND MOLECULAR LIFE SCIENCES*, **61**(2004):1198-207.
4. Buliard A, Fortenfant F, Ghillani-Dalbin P, Musset L, Oksman F, Olsson NO, Analysis of nine autoantibodies associated with systemic autoimmune diseases using the Luminex technology. Results of a multicenter study. *Ann Biol Clin (Paris)*. , **63**(2005):51-8.
5. Carson RT, Vignali DAA, Simultaneous quantitation of 15 cytokines using a multiplexed flow cytometric assay. *Journal of Immunological Methods*, **227** (1999): 41-52.
6. Hindson BJ, McBride MT, Makarewicz AJ, Henderer BD, Setlur US, Smith SM, Gutierrez DM, Metz TR, Nasarabadi SL, Venkateswaran KS, Farrow SW, Colston BW Jr, Dzenitis JM, Autonomous detection of aerosolized biological agents by multiplexed immunoassay with polymerase chain reaction confirmation. *Anal Chem*, **77**(2005):284-9.
7. Wilson WJ, Erler AM, Nasarabadi SL, Skowronski EW, Imbro PM, A multiplexed PCR-coupled liquid bead array for the simultaneous detection of four biothreat agents. *Molecular and Cellular Probes*, **19** (2005): 137-144.
8. Biagini RE, Sammons DL, Smith JP, MacKenzie BA, Striley CA, Robertson SA,

Snawder JE, Quinn CP, 2005, Simultaneous measurement of specific serum IgG responses to five select agents. *Anal Bioanal Chem.*, 382(2005):1027-34

9. Barad O, Meiri E, Avniel A, Aharonov R, Barzilai A, Bentwich I, Einav U, Gilad S, Hurban P, Karov Y, Lobenhofer EK, Sharon E, Shibolet Y, Shtutman M, Bentwich Z, Einat P, MicroRNA expression detected by oligonucleotide microarrays: system establishment and expression profiling in human tissues. *GENOME RESEARCH*, 14(2004):2486-94.

10. Hadd AG, Laosinchai-Wolf W, Novak CR, Badgett MR, Isgur LA, Goldrick M, Walkerpeach CR, Microsphere Bead Arrays and Sequence Validation of 5/7/9T Genotypes for Multiplex Screening of Cystic Fibrosis Polymorphisms. *Journal of Molecular Diagnostics*, 6(2004): 348-355.

11. Hurley JD, Engle LJ, Davis JT, Welsh AM, Landers JE, A simple, bead-based approach for multi-SNP molecular haplotyping. *Nucleic Acids Research*, 32 (2005): e186.

12. Dias D, Van Doren J, Schlottmann S, Kelly S, Puchalski D, Ruiz W, Boerckel P, Kessler J, Antonello JM, Green T, Brown M, Smith J, Chirmule N, Barr E, Jansen KU, Esser MT, 2005, Optimization and validation of a multiplexed luminex assay to quantify antibodies to neutralizing epitopes on human papillomaviruses 6, 11, 16, and 18. *Clin Diagn Lab Immunol*, 12 (2005):959-69.

13. Johnson AJ, Noga AJ, Kosoy O, Lanciotti RS, Johnson AA, Biggerstaff BJ, Duplex microsphere-based immunoassay for detection of anti-west nile virus and anti-st. Louis encephalitis virus immunoglobulin m antibodies. *Clin Diagn Lab Immunol*, 12 (2005):566-74.

14. Khan IH, Kendall LV, Ziman M, Wong S, Mendoza S, Fahey J, Griffey SM, Barthold SW, Luciw PA., Simultaneous serodetection of 10 highly prevalent mouse infectious pathogens in a single reaction by multiplex analysis. *Clinical and Diagnostic Laboratory Immunology*, 12 (2005): 513-519.

15. Lukacs Z, Dietrich A, Ganschow R, Kohlschutter A, Kruithof R, Simultaneous determination of HIV antibodies, hepatitis C antibodies, and hepatitis B antigens in dried blood spots--a feasibility study using a multi-analyte immunoassay. *Clin Chem*

Lab Med, 43 (2005):141-5.

16. Schmitt M, Bravo IG, Snijders PJ, Gissmann L, Pawlita M, Waterboer T, Bead-based multiplex genotyping of human papillomaviruses. *J Clin Microbiol.*, 44(2006):504-12.

17. Lash GE, Scaife PJ, Innes BA, Otun HA, Robson SC, Searle RF, Bulmer JN, Comparison of three multiplex cytokine analysis systems: Luminex, SearchLighttrade mark and FAST Quant(R). *J Immunol Methods.*, 309(2006):205-8.


18. Amos JA, Bridge-Cook P, Ponek V, Jarvis MR, A universal array-based multiplexed test for cystic fibrosis carrier screening. *Expert Rev Mol Diagn.*, 6(2006):15-22.

19. Das S, Brown TM, Kellar KL, Holloway BP, Morrison CJ, DNA probes for the rapid identification of medically important *Candida* species using a multianalyte profiling system. *FEMS Immunol Med Microbiol*, 46(2006):244-50.

20. Prince HE, Lape-Nixon M, Matud J, Evaluation of a Tetraplex Microsphere Assay for *Bordetella pertussis* Antibodies. *Clin Vaccine Immunol.*, 13(2006):266-70.

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