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Abstract

Shigella sonnei, a primary travel-related pathogen of diarrheal diseases in Taiwan, occasionally causes large scale disease outbreak in crowded institutes. It is also one of the most common foodborne pathogens in North America. Since *Shigella* is transmitted from an infected person to another usually by a fecal-oral route, it causes child mortality due to diarrhea in countries or regions with poor sanitary condition.

Pulsed-field gel electrophoresis (PFGE) is the standard molecular subtyping tool for PulseNet, a molecular subtyping network for foodborne disease surveillance. It is highly discriminatory but, in times, it is insufficient in discriminatory some epidemiologically unrelated isolates of monomorphic bacterial species. PFGE is also costly, labor-intensive and time-consuming and generates image data that is disadvantageous to pattern comparison among laboratories. Multilocus variable number tandem repeat analysis (MLVA) consisting only a small set of VNTR loci can have high discriminatory power over PFGE. MLVA generates simple integer data that are easy to analyze and compare among laboratories. However, to be a standard subtyping tool for disease surveillance, a MLVA protocol needs to be validated. In this study, an MLVA protocol is validated by laboratories of Japan, Korea, Malaysia and India. The phase I validation decided a common MLVA protocol consisting of 7 VNTR loci (SS1, SS3, SS6, SS9, SS10, SS11, SS13). The 7 VNTR loci can generate at least 8.5×10^7 genotypes that provide sufficient discriminatory power for subtyping *S. sonnei* isolates for disease surveillance and outbreak investigation.

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副本: ohnishi7@nih.go.jp

Dear Dr. Chiou,

We have now an possible outbreak of *Shigella sonnei*.
But we have not yet identified the source of contamination.

I would very much appreciate if you examine our MLVA data for your database, and let me know when you find some isolates of match, or single- or double-locus variants in your database;

SS1 9

SS10v 3

SS11v 4

SS13v 4

SS3 15

SS6 13

SS9v 11

The values indicate the repeat number for the loci.

Best regards,

Hidemasa Izumiya

--

Hidemasa Izumiya, Ph.D.

Chief, Laboratory of Unit II (Enteric Infection),
Department of Bacteriology,
National Institute of Infectious Diseases

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附錄一、Laboratory Standard Operating Procedure for MLVA of *Shigella sonnei* –
Applied Biosystems Genetic Analyzer 3130 XL Platform

1. PURPOSE: to describe the standardized laboratory protocol for molecular subtyping of *S. sonnei*
2. SCOPE: to provide the PulseNet participants with a single protocol for performing MLVA of *S. sonnei*, thus ensuring inter-laboratory comparability of the generated results.
3. DEFINITIONS:
 - 3.1 MLVA: Multilocus variable-number tandem repeat analysis
 - 3.2 VNTR: Variable-number tandem repeat
 - 3.3 DNA: Deoxyribonucleic acid
 - 3.4 DNase: Deoxyribonuclenase
 - 3.5 PCR: Polymerase chain reaction
 - 3.6 dNTP: Deoxyribonucleotide triphosphate
 - 3.7 SOP: Standard Operating Procedure
4. RESPONSIBILITIES/PROCEDURE
 - 4.1. Biosafety warning: *S. sonnei* is a human pathogen capable of causing serious disease. Always use a minimum of Biosafety Level 2 practices and extreme caution when transferring and handling strains of this pathogen. Work in a biological safety cabinet when handling large amounts of cells. Disinfect or dispose of all plastic ware and glassware that come in contact with the cultures in a safe manner.
 - 4.2. Reagents, supplies and equipment needed for DNA template preparation
 - 4.2.1 Trypticase soy agar (TSA) or comparable media
 - 4.2.2 Inoculation Needle (Copan, Code# 176CSR40)
 - 4.2.3 1.5 ml microcentrifuge tubes
 - 4.2.4 DNase-free, molecular biology -grade water
 - 4.2.5 Vortex
 - 4.2.6 Boiling water bath or Dry bath accommodating 1.5 ml tubes
 - 4.2.7 Tabletop centrifuge
 - 4.2.8 Pipets (200 µl) for aliquoting 100 µl of DNase-free, molecular biology-grade water

4.2.9 Filtered sterile pipet tips

4.3. Reagents, supplies and equipment needed for PCR

4.3.1 DNA lysates from isolates (keep at -20°C or -80°C freezer for long term)

4.3.2 PCR primers (see Appendix Table 1)

4.3.2.1 Fluorescent-labeled forward primers

4.3.2.2 Unlabeled reverse primers

4.3.2.3 Applied Biosystems sells primers labeled with the four dyes needed for the protocol

4.3.2.4 Divide the concentrated stocks in portions and store at -80°C freezer

4.3.2.4.1 One vial should contain enough to prepare 25-50 µl of working solution. Avoid repetitive freeze-thaw cycles of concentrated primer stocks.

4.3.2.5 The working solutions can be stored at either -20°C or -80°C freezer

4.3.2.6 Prepare primer mix solutions (100 µl) as follows:

Primer mix (100 µl)	Locus	Primer/H2O	Working solution (µM)	Volume (µl)	Final conc. in mix (µM)	Final conc. in PCR (µM)
M1	SS1	SS1-F	100	1	1	0.05
		SS1-R	100	1	1	0.05
	SS3	SS3-F	100	2	2	0.1
		SS3-R	100	2	2	0.1
	SS6	SS6-F	100	4	4	0.2
		SS6-R	100	4	4	0.2
	SS9	SS9-F	100	2	2	0.1
		SS9-R	100	2	2	0.1
			H2O		82	
M2	SS10	SS10-F	100	1	1	0.05
		SS10-R	100	1	1	0.05
	SS11	SS11-F	100	1	1	0.05
		SS11-R	100	1	1	0.05
	SS12	SS12-F	100	5	5	0.25
		SS12-R	100	5	5	0.25

SS13	SS13-F	100	1	1	0.05
	SS13-R	100	1	1	0.05
	H2O		84		

4.3.2.7 Divide primer mixes into aliquots and store at either -20°C or -80°C freezer. Use new primer mixes every month or if a significant drop in the fluorescence level is observed

4.3.3 96-well polypropylene PCR plates (Axygen, Cat. No. PCR-96M2-HS-C) or Microamp PCR tubes without caps (Applied Biosystems, Cat. No. N8010533)

4.3.4 8-well strip caps for the polypropylene plate Axygen, Cat. No. PCR-02CP-C) or MicroAmp strip caps for the individual tubes (Applied Biosystems, Cat. No. N8010535)

4.3.5 DNase-free, molecular biology-grade water

4.3.6 1.5 ml Eppendorf microcentrifuge tubes

4.3.7 dNTP mixture (dNTP mix containing all four nucleotides; BioVan, Cat. No. BV105)

4.3.8 Taq Polymerase and 10X buffer (BioVan, Cat. No. TQ1000)

4.3.9 25 mM MgCl₂ (Roche, Cat. No. 11699113001)

4.3.10 PCR Cooling block (VWR International, Cat. No. 62111-762)

4.3.11 DNA Engine (Biorad), GeneAmp (Applied Biosystems) or similar thermocycler with a heated lid option and 96-well block format

4.3.12 Parafilm M, 4" width (VWR, Cat. No. 52858-032)

4.3.13 Complete set (1000 µl, 200 µl, 100 µl, 20 µl, 10 µl and 2 µl) of single channel pipettors for mastermix set-up ("clean set")

4.3.14 10 µl and 100 µl single channel pipettors for preparing and adding DNA templates

4.3.15 Filtered tips for pipettors

4.4. Reagents, supplies and equipment needed for Genetic Analyzer 3130XL

4.4.1 DNase-free, molecular biology-grade water

4.4.2 PCR Cooling block (VWR International, Cat. No. 62111-762)

4.4.3 10 µl, 100 µl, and 1000 µl single channel pipettors

4.4.4 10 µl and 200 µl multichannel pipettors

4.4.5 Filtered pipette tips

4.4.6 Sterile solution basins

4.4.7 1.5 ml Eppendorf microcentrifuge tubes

4.4.8 96-well polypropylene PCR plates (for dilutions; Axygen, Cat. No.

PCR-96M2-HS-C)

- 4.4.9 MicroAmp Optical 96-well reaction plates (Applied Biosystems, Cat. No. N801-0560)
 - 4.4.10 96-well plate base (Applied Biosystems, Cat. No. 4317237)
 - 4.4.11 Rubber septa for 96-well reaction plates (Applied Biosystems, Cat. No. 4315933)
 - 4.4.12 96-well plate retainer (Applied Biosystems, Cat. No. 4317241)
 - 4.4.13 Hi-Di Formamide, 25 ml (Applied Biosystems, Cat. No. 4311320)
 - 4.4.14 *GeneScan 500LIZ Size Standard* (Applied Biosystems, Cat. No. Number 4322682)
 - 4.4.15 Multi-Capillary DS-33 (Dye Set G5) Matrix Standard Kit (Applied Biosystems, Cat No. 4345833)
 - 4.4.15.1 Needed periodically to establish the system dye color spectra for the instrument
 - 4.4.16 3130XL & 3100 Capillary Array 36 cm (Applied Biosystems, Cat. No. 4333464)
 - 4.4.17 10X Running Buffer with EDTA, 25 ml (Applied Biosystems, Cat. No. 402824)
 - 4.4.18 3130XL POP7 Performance Optimized Polymer, 3.5 ml (Applied Biosystems, Cat.No. 4363785)
 - 4.4.19 Parafilm M, 4" width (VWR, Cat. No. 52858-032)
 - 4.4.20 A rinse bottle containing distilled water
 - 4.4.21 Centrifuge with a microtiter plate rotor
 - 4.4.22 Thermal cycler or heating block for denaturation
- 4.5. Isolate and DNA template preparation
- 4.5.1 Day 0:
 - 4.5.1.1 Streak an isolated colony from test cultures to TSA plate (or comparable media). Incubate cultures at 37°C for 14-18 hrs.
 - 4.5.2 Day 1:
 - 4.5.2.1 For each isolate to be typed, aliquot 100 µl of sterile, molecular biology-grade water into 1.5 ml microfuge tubes. Use a sterile, disposable inoculation needle to scratch 1 colony or 1 µl loop to pick 2-3 colonies (about half of a loopful); stir the inoculation needle or agitate/rotate the loop in a microfuge tube to release the bacteria into the water. Cap and vortex for 10-15 seconds to disperse any clumps.
 - 4.5.2.2 Place the tubes in a 100°C heat block or water bath for 10

minutes. Cool briefly on ice or in fridge and centrifuge for 1 minutes at 15,000 rpm. Place on ice or in fridge while preparing PCR reactions. These DNA templates can be stored at - 20°C or -80°C for several years.

4.6. PCR procedure

4.6.1 Day 1:

4.6.1.1 Fill out, save with the run name, and print a PCR template worksheet (see Appendix Table 2) with appropriately labeled samples (a maximum of 46 unknown isolates/plate + one positive control [DNA of strain C07.0755] and one negative control)

4.6.1.1.1 For each isolate, two wells must be labeled as follows: “ID_M1” where “ID_M1” represents the isolate ID (be sure to use the exactly same isolate ID that is used in the PFGE gels uploaded to BioNumerics database) and “M1” represents one of the two specific multiplex PCR reactions (M1, M2).

4.6.1.2 Fill out, and print a PCR mastermix calculation worksheet (see Appendix Table 3) by typing the number of isolates to be tested (plus 2-3 extra) in the PCR mastermix calculators labeled M1 and M2. This number is highlighted in RED and is next to “number of samples to be analyzed”. The mastermixes for reactions M1 and M2 for one sample are as follows:

Reagent	Volume (µl)	Final conc.
PCR water	4.9	—
DMSO (10%)	1	1%
PCR buffer (10 X)	1	1X
MgCl ₂ (25 mM)	0.6	3 mM*
dNTP (2.5 mM)	0.8	0.2 mM
Primer mix (M1 or M2)	0.5	various
Tag (5 u/µl)	0.2	1 U/9µl
Total	9	

*1X Buffer contains 1.5 mM MgCl₂

4.6.1.2.1. NOTE: these primer concentrations serve as a starting point. Since laboratory-specific factors, such as the age of the primer stocks, calibration status of the thermocyclers

and pipettes, etc. affect amplification efficiency, each laboratory will have to re-optimize the primer concentrations for optimal detection of all targets. However, any other parameters stated in the SOP should not be changed.

4.6.1.3 Thaw all reagents and supplies needed for PCR reactions and place on ice; keep primers light protected as much as possible

4.6.1.3.1 NOTE: PCR mastermixes should be set up in a clean hood that is dedicated just for this purpose and where no cultures or DNA are handled

4.6.1.4 Prepare the two separate PCR mastermixes in 1.5 ml Eppendorf tubes following the instructions in the PCR mastermix calculation worksheet (see Appendix Table 3). Keep the mastermix on ice while preparing. Add the mastermix components in the following order: water, DMSO, 10x PCR buffer, MgCl₂, dNTPs, primer mix (M1 or M2), and then finally Taq polymerase. Mix the reaction mixture by vortexing briefly.

4.6.1.4.1 NOTE: All components except Taq polymerase should be vortexed thoroughly before adding to the mastermix. Taq may be briefly centrifuged, if necessary to pull the enzyme down to the bottom of the tube

4.6.1.5 Place a 96-well PCR plate or required number of PCR tubes in a PCR cooling block

4.6.1.6 Dispense 9.0 µl of each mastermix into the appropriate rows of the 96-well polypropylene plate / PCR tubes as noted in the PCR template worksheet (see Appendix Table 2)

4.6.1.7 Add 1.0 µl of DNA lysate to each of the two different wells representing the two PCR reactions. Use *S. sonnei* strain C07.0755 as a positive control and water as a negative control.

4.6.1.8 Cover all wells / tubes with 8-well strip caps and firmly clamp down to avoid any evaporation during PCR amplification

4.6.1.9 Recommendation: briefly spin down the plate / tubes to remove any air bubbles

4.6.1.10 Setup the following PCR parameters and save the program as "SSMLVA"

* 94°C for 5 min Step 1

- * 94°C for 45 sec Step 2
- * 55°C for 50 sec Step 3
- * 72°C for 60 sec Step 4
- * Go to step 2, 29x Step 5
- * 72°C for 5 min Step 6
- * Indefinite hold at 25°C Step 7

4.6.1.10.1 NOTE: Make sure to use the heated lid option on the PCR block and tube (calculated) temperature control

4.6.1.11 When the PCR is complete store the amplification products light-protected at 4°C until ready to run on the sequencer. If the fragment analysis is not performed the same day, the plate should be stored at -20°C or -80°C. The PCR products are stable for approximately one month, when stored frozen.

4.7. Initial setup of Genetic Analyzer 3130XL instrument:

4.7.1 Turn on the computer, which is connecting to 3130XL, set up Gene Mapper v4.0

4.7.1.1 Double-click on the shortcut icon for GeneMapper v.4.0 and enter the appropriate password to access the software. The main menu window will open

4.7.1.2 Set up the size standard:

4.7.1.2.1 From the “Tools” drop-down menu, select “GeneMapper Manager”. A “GeneMapper Manager” window will open.

4.7.1.2.2 Select the “Size Standards” tab

4.7.1.2.3 Click on “New”, leave the default option “Basic or Advanced” checked and click “OK”. “Size Standard Editor” window will appear.

4.7.1.2.4 Name the new size standard “500LIZ”

4.7.1.2.5 Leave the default option “Orange” as “Size Standard Dye”

4.7.1.2.6 Enter sizes for each peak in the table

4.7.1.2.7 When finished, click “OK”

4.7.1.3 Set up analysis method:

4.7.1.3.1 From the “Tools” drop-down menu, select “GeneMapper Manager”. A “GeneMapper Manager” window will open.

4.7.1.3.2 Select the “Analysis Methods” tab

4.7.1.3.3 Click on “New”, leave the default option “Microsatellite” checked and click “OK”. “Analysis Methods Editor” window

will appear.

4.7.1.3.4 Name the new method “SSMLVA” and click “OK”

4.7.1.3.5 Highlight the new method name “SSMLVA” and click “Open”

4.7.1.3.6 Select the “Peak Detector” tab and change “Peak Detection Algorithm” to “Advanced” from the drop-down menu

4.7.1.3.7 Input the following analysis settings and click “OK” when finished:

4.7.1.3.7.1 Analysis: Full Range

4.7.1.3.7.2 Sizing: Partial Sizing

4.7.1.3.7.2.1 Start Size: 60

4.7.1.3.7.2.2 Stop Size 600

4.7.1.3.7.3 Smoothing: none

4.7.1.3.7.4 Baseline window: 51 pts

4.7.1.3.7.5 Size Calling Method: Local Southern Method

4.7.1.3.7.6 Peak Amplitude Thresholds

4.7.1.3.7.6.1 B: 100

4.7.1.3.7.6.2 G: 100

4.7.1.3.7.6.3 Y: 50

4.7.1.3.7.6.4 R: 100

4.7.1.3.7.6.5 O: 50

4.7.1.3.7.7 Min. Peak Half Width: 2 pts

4.7.1.3.7.8 Polynomial Degree: 3

4.7.1.3.7.9 Peak Window Size: 15 pts

4.7.1.3.7.10 Slope Threshold

4.7.1.3.7.10.1 Peak Start: 0.0

4.7.1.3.7.10.2 Peak End: 0.0

4.7.2 Turn on the power of Analyzer, wait until the Status Light turns to Green. Click on the “Run 3130XL Data Collection v 3.0” icon of the computer. The “Service Console” window will appear. If the connections are functioning properly, the “Messaging Service”, “Data Service”, “Instrument Service” and “Viewer” icons will change from red circles to green squares. The main window of the “Foundation Data Collection” software will open. Check to make sure system status is green.

4.7.3 NOTE: steps 4.7.4 and 4.7.5 only need to be performed before the very first run

4.7.4 Set up a results group:

4.7.4.1 On the left side of the window under “GA Instruments”, highlight “Results group” and click on “New”.

4.7.4.2 The “Results Group Editor” window will appear

4.7.4.3 Under “General” tab, name the new results group as “SSresults”

4.7.4.4 Under “Analysis” tab, select “GeneMapper-Generic” from the “Analysis type” drop-down menu. Leave “Analysis actions” unchecked

4.7.4.5 Under “Destination” tab, keep the default root destination (E:\AppliedBiosystems\udc\DataCollection\Data)

4.7.4.6 Under “Naming” Tab, select “Sample Name” from the first “Sample File NameFormat” drop-down menu, and select “Plate Name” from the first “Run Folder Name Format” drop-down menu. Leave all other fields blank.

4.7.4.7 Under “Automated Processing” tab, leave the default “Only when the results group is complete” checked

4.7.5 Set up running conditions:

4.7.5.1 Under “GA Instruments”, click on + to expand “ga3130XL” subfolders

4.7.5.2 Highlight “Protocol Manager”. The instrument protocols will be listed on the right side of the window.

4.7.5.3 Click on “New”, and the “Protocol Editor” window will appear

4.7.5.3.1 Name the new protocol “SSTest”

4.7.5.3.2 Leave the protocol type as “Regular”

4.7.5.3.3 From the “Run Module” drop-down menu, select “Fragment Analysis 36_POP7_1” as running conditions

4.7.5.3.3.1 NOTE: These are the instrument default running conditions for 36 cm capillary array and POP7 polymer. You can check the running conditions and modify them by highlighting “Module Manager” subfolder under ga3130XL and by double clicking on the protocol name. The default conditions are:

4.7.5.3.3.1.1 Oven_Temperature: 60°C

4.7.5.3.3.1.2 Poly Fill Vol: 6500 Steps

4.7.5.3.3.1.3 Current Stability: 5.0 μ Amps

4.7.5.3.3.1.4 PreRun_Voltage: 15.0 kVolts

4.7.5.3.3.1.5 Pre-Run-Time: 180 sec.

- 4.7.5.3.3.1.6 Injection_Voltage: 1.2 kVolts
- 4.7.5.3.3.1.7 Injection_Time: 23 sec.
- 4.7.5.3.3.1.8 Voltage_Number_Of_Steps: 20 nk
- 4.7.5.3.3.1.9 Voltage_Step_Interval: 15 sec.
- 4.7.5.3.3.1.10 Data_Delay_Time: 60 sec.
- 4.7.5.3.3.1.11 Run_Voltage: 15.0 kVolts
- 4.7.5.3.3.1.12 Run time: 1200sec.

4.7.5.3.4 From the “Dye Set” drop-down menu, select “G5” as dye set

4.8. Genetic Analyzer 3130XL instrument preparation before each run

4.8.1. Day 1

4.8.1.1 Make sure service console is open and all components are green. Make sure a capillary array is installed in the instrument. For installation, follow the instructions of the “Install Array Wizard”. You can find the wizards by expanding the 3130XL subfolder and by highlighting “Manual Control”. The “Wizards” drop-down menu will appear on top of the window.

4.8.1.2 Setup the plate run on the 3130XL

4.8.1.2.1 Under the ga3130XL, highlight “Plate Manager”

4.8.1.2.2 Click on “New” and a “New Plate Dialog Window” will appear

4.8.1.2.3 Name the run following a standardized naming system: use the unique identifier code that was assigned to operator of laboratory for the first two letters of the file name. The next two spaces will indicate the year and the next four spaces will indicate the month and the date the run was performed. For example YW100426 is a run made by an operator “YW” on April 26th 2010. If several runs are performed the same day, separate the file names by using sequential numbers, for example YW100426-1, YW100426-2.

4.8.1.2.4 From the “Application” drop-down menu, select “GeneMapper-Generic”

4.8.1.2.5 Type in the Owner Name and the Operator Name

4.8.1.2.6 Click “OK” and the “GeneMapper Plate Editor” window will appear

4.8.1.2.6.1 Type in the sample IDs

4.8.1.2.6.2 Select "SSresults" from the "Results Group 1" drop-down menu for the first sample, highlight the results group column, and select "Fill Down" from the "Edit" drop-down menu

4.8.1.2.6.3 Select "SSTest" from the "Instrument Protocol 1" drop-down menu, highlight the instrument protocol column, and select "Fill Down" from the "Edit" drop-down menu

4.8.1.2.6.4 Select "500LIZ" from "Size standard"

4.8.1.2.6.5 Select "SSMLVA" from the Analysis Method

4.8.1.2.6.6 Select "None" from Panel

4.8.1.2.6.7 Select "Sample" from Sample Type

4.8.1.2.6.8 Click "OK"

4.8.1.3 Install the POP7 polymer in the instrument

4.8.1.3.1 Expand the 3130XL subfolder and highlight "Manual Control". The "Wizards" drop-down menu will appear on top of the window.

4.8.1.3.2 Follow the instructions of the "Replenish Polymer Wizard"

4.8.1.3.2.1 NOTE1: if an old polymer (been in the instrument > 7 days) is switched to a new one follow the instructions of the "Water Wash Wizard" instead. Use distilled water that has been heated to 37-40°C for water wash.

4.8.1.3.2.2 NOTE2: water wash should be performed once a week

4.9. Fragment analysis sample preparation

4.9.1 Day 1

4.9.1.1 Thaw the Hi-Di Formamide, the 500LIZ DNA size standard and place on ice.

4.9.1.1.1 NOTE: aliquot Hi-Di Formamide (500 µl/tube) and the size standard (50 µl/tube) in order to avoid frequent freeze-thaw cycles.

4.9.1.2 Prepare a 96-well plate for diluting the PCR reactions. Using a 200 µl multichannel pipettor and a solution basin, dispense 19 µl of molecular-grade water in the required number of wells.

4.9.1.3 Remove the plate / tubes with the PCR reactions from the thermocycler. Briefly spin down the plate / tubes, if necessary.

Use a 10 µl multichannel pipettor to transfer 1 µl of each PCR reaction directly across to the corresponding set of wells in the dilution plate. In order to avoid cross-contamination, remove the strip cap from just one column at a time and recap the column before opening the next one.

4.9.1.4 Using a 200 µl multichannel pipettor, mix the dilutions by pipetting up and down a few times. Cover the plate with parafilm and put in the fridge or on ice.

4.9.1.5 Prepare a fragment analysis master mix containing DNA size standard and Hi-Di Formamide for the samples. The fragment analysis mastermix can be calculated as follows. Vortex briefly and place on ice.

Reagents	Frag. anal. mastermix
Hi-Di Formamide	9 µl x (# samples +3)=
500LIZ size standard	0.1 µl x (# samples +3)=

4.9.1.6 Place a MicroAmp Optical 96-well sample plate in a cold block. Aliquot 9 µl of the prepared fragment analysis mastermix to the required number of wells. Cover the plate loosely with Parafilm.

4.9.1.7 Using the 10 µl multichannel pipettor, add 1 µl of 1:20 diluted PCR reactions to the appropriate columns in the sample plate. Keep sliding the Parafilm sheet from column to column to keep track of the sample order.

4.9.1.8 Denature templates by heating the reaction plate in a thermocycler at 95°C for 3 min, and place on ice immediate.

4.9.1.9 While the templates denatured, prepare 30 ml 1X running buffer with EDTA and turn the oven of the 3130XL on using Manual Control. Change set point to 60°C.

4.9.1.9.1 Add 3 ml of 10X Running Buffer with EDTA in a 50 ml graduated cylinder

4.9.1.9.2 Add purified water to bring the total volume up to 30 ml

4.9.1.9.3 Mix well by inverting

4.9.1.10 Briefly spin down the sample plate to remove any air bubbles

4.9.1.11 Place the sample plate in a plate base. Seal the plate with the rubber septa. Snap the plate retainer onto the plate and the plate base.

4.9.1.12 Place the plate assembly and the buffer into the Genetic Analyzer 3130XL.

4.9.1.12.1 Push on the “Tray” button at the front of the Genetic Analyzer 3130XL to bring the autosampler to the forward position. Open the instrument doors.

4.9.1.12.2 Place the plate assembly on the autosampler in position A or B with the notched end of the plate base away from you

4.9.1.12.3 Add 1X running buffer to the anode and cathode reservoirs

4.9.1.12.4 Add distilled water to the waste and rinse reservoirs

4.9.1.12.5 Close the instrument doors and wait for the green light to illuminate 4.9.1.13. Start the run after making sure the system status is green.

4.9.1.13.1 Expand 3130XL computer name icon.

4.9.1.13.2 Highlight “Run Scheduler”

4.9.1.13.3 Click on “Find All”

4.9.1.13.4 Highlight the plate name for the run

4.9.1.13.5 Link the plate by clicking on the yellow squares on the right side of the window

4.9.1.13.6 The “Start Run” button will turn green indicating that the run can be started. Click on this button, and then click “OK” on the alert window. The run will start.

4.10. Viewing and exporting data from the Genetic Analyzer 3130XL

4.10.1 Day 2

4.10.1.1 Double-click on the shortcut icon for GeneMapper v.4.0 and enter the appropriate password to access the software. The main menu window will open.

4.10.1.2 From the “File” drop-down menu select “Add Samples to Project”

4.10.1.3 Find the folder containing the data file to be analyzed: My Computer → E: → Applied Biosystems → UDC → DataCollection → Data

4.10.1.4 Highlight the desired file(s) and click on “Add to List”. File(s) will appear in the window on the right. Click “Add & Analyze” below the file list to return to the original screen

4.10.1.5 Name the project with the run name (for example, “YW100426”)

and click “OK”

4.10.1.6 A successful analysis is indicated by green squares. Yellow triangles indicate problematic components (i.e. missing size standard peaks). Red circles indicate results fell below acceptable quality values. Samples with yellow or red circles in the SQ column should be selected for re-analysis.

4.10.1.6.1 To resolve failed analyses due to sub-optimal molecular marker peak profile (i.e. miscalling of peaks), select a row with a yellow triangle or red circle in the SQ column and click on the “Size Match Editor” icon on the toolbar. “Size Match Editor” view will appear.

4.10.1.6.2 Place the cursor near the X-axis to activate the magnifying lens, and then pull up (mouse left-click and hold) to zoom in a specific area to facilitate editing.

4.10.1.6.3 Left-click at the base of a peak to select. Right-click and select “Add”, “Delete”, or “Change”.

4.10.1.6.4 Select the correct molecular weight for the selected peak from the drop-down menu. Repeat this process for all other miscalled peaks. Click “OK” when finished.

4.10.1.6.5 After the size standard has been adjusted, click the “Play” button to re-analyze the data. After a successful analysis, the samples will have green squares under the SQ column. If the size standard cannot be adjusted, the reaction is considered a fragment analysis failure and must be re-run.

4.10.1.7 Check the fragment result data (the fluorescent peaks) for each well by highlighting the well ID and by clicking on the “Display Plots” icon on the toolbar. Make sure that all eight VNTRs amplified in the C07.0755 control and that the fragment size is within ± 1 bp from that specified in the Appendix Table 4. Write down any failed reactions in the PCR/ABI worksheet. Make a note of non-specific bands and primer-dimers.

4.10.1.8 Export the peak file:

4.10.1.15.2 NOTE: the following columns should appear in the exported table in the following order from left to right: “Dye/Sample”, “Sample file name”, “Marker”, “Size”, “Height”, “Area”, “Data point”. The format of the table can

be modified by selecting “Table setting editor” from the “Tools” drop down menu. Select the “Genotype” tab and make sure that the boxes for the above mentioned columns are checked and no additional boxes are checked.

4.10.1.8.2 Highlight the samples for which you want to export peak data

4.10.1.8.3 Click on “Display Plots” icon on the toolbar

4.10.1.8.4 Click on “Sizing Table” icon on the toolbar and a table will appear below the electropherograms

4.10.1.8.5 From the “File” drop-down menu, select “Export Table”

4.10.1.8.6 Select the location (for example a flash drive) where you want to export the data

4.10.1.8.7 Name the export file with the run name (for example YW100426) and make sure the file type is tab-delimited text (.txt) file

4.10.1.9 The remaining gel can stay in the instrument if it is going to be used within seven days

4.10.1.10 Convert the sizes of fragments into copy numbers (see Excel file: *S. sonnei* worksheet_Copy number and amplicon sizes)

4.10.1.11 Note: Sequencing is recommended to confirm the copy number of an amplicon with an unexpected size. The flanking regions for the VNTRs used in this protocol are listed in the Excel file: *S. sonnei* worksheet_Flanking regions

5. RREFERENCES

- 5.1 Liang SY, Watanabe H, Terajima J, Li CC, Liao JC, Tung SK, Chiou CS: **Multilocus Variable-Number Tandem Repeat Analysis for Molecular Typing of *S. sonnei***. *Journal Clinical Microbiology* 2007, **45**(11):3574-3580

6. CONTACTS

6.1 Chien-Shun Chiou, Ph.D. Central Region Laboratory, Center for Research and Diagnostics, Centers for Disease Control, Taiwan. TEL: 886-4-24750452; Email: nipmcsc@cdc.gov.tw

7. AMENDMENTS:

Appendix Table 1.

S. sonnei MLVA PCR primers.

Locus	Dye_Forward primer (5' to 3')	Reverse primer (5' to 3')
M1		
	SS1 VIC_TTGCCAGTACACCTCACTCG	GCGTCGGCGTTAATATCACT
	SS3 6-FAM_CTGGGAGATGAACAGGAGGA	ATGCCAGCGACAAGTTTCTT
	SS6 NED_GAGTCGCTAAACGCTTGCTT	GGGAAATAGAGCGGACCTTT
	SS9 PET_CGCAATCAGCAAAACAAAGA	GCGATGCTGGAAAACTGAT
M2		
	SS10 6FAM_ACGGTGGGCTTTCTCTACCT	GCGAGGGGGAGATCAGTATT
	SS11 VIC_CTGGTCCGGGAGATTATCG	CTGTTTCAGCGGTCTCTTCC
	SS12 NED_GCTGTAGGCACGGAAAAGAA	TGGATATTGTGCAGGGTTCA
	SS13 PET_AGACGCTGGCTTATGACGAT	GCCACTGTCCTGAATGGTCT

Appendix Table 2

S. sonnei MLVA PCR Template and ABI worksheet (see Excel file: *S. sonnei* worksheet_PCR template)

PCR run date:

PCR instrument:

ABI run date:

ABI instrument:

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
J												

Appendix Table 3

Mastermixes for PCR M1 and M2 (see Excel file: *S. sonnei* worksheet_Mastermix calculator)

Technician:

<i>S. sonnei</i> multiplexes PCR M1 & M2			
	20	Number of (X) samples to be analyzed	
	10	µl Total Reaction Volume	
	Volume (µl) for 1 sample	Conc. in PCR reaction	Volume (µl) for X samples
PCR water	4.9	—	98
DMSO (10%)	1	1%	20
PCR buffer (10 X)	1	1X	20
MgCl ₂ (25 mM)	0.6	1.5 mM	12
dNTP (2.5 mM)	0.8	0.2 mM	16
Primer mix M1 or M2	0.5		10
Tag (5 u/µl)	0.2	0.1 u/µl	4
Total	9		180

Appendix Table 4

Sizes of VNTRs for positive control strain C07.0755

Multiplex	Locus	Repeat_bp	Copy no.	Size_amplicon_bp	Size estimated by		
					500LIZ*	600LIZ	ROX625
M1	SS1	7	15	230	216.01	213.95	216.48
	SS3	7	17	233	230.02	227.65	231.29
	SS6	7	24	320	314.83	311.44	317.5
	SS9	6	9	215	215.6	213.31	216.11
M2	SS10	6	3	198	194.61	191.6	194.39
	SS11	6	4	199	197.21	194.33	197.04
	SS12	9	2	188	184.93	181.54	184.74
	SS13	6	3	236	236.73	234.79	238.73

* Marker 250 bp is not included as recommended by the Supplier.